



Pancreatic Insulin Content by Acid-Ethanol Extraction

Version: 1

Replaced by version: N/A

Edited by: Ed Leiter, Ph.D. Jackson Laboratories

Summary: Used to calculate the insulin content in the pancreas.

Protocol:

- 1) $\frac{1}{4}$ - $\frac{1}{2}$ of the pancreas is placed into 5 ml Acid-Ethanol (1.5% HCl in 70% EtOH) in a 15 ml conical vial.
- 2) Incubate O/N at -20°C .
- 3) Homogenize tissue (I use a Polytron homogenizer).
- 4) Incubate O/N at -20°C .
- 5) Centrifuge at 2000 rpm 15 min at 4°C (Sorvall RT6000).
- 6) Transfer aqueous solution to a new 15 ml conical vial.
- 7) Neutralize 100 μl of Acid-Ethanol extract with 100 μl 1M Tris pH 7.5.
- 8) Dilute further (1:100, 1:1000, or 1:5000 depending upon the strain) in Insulin ELISA sample diluent.
- 9) Run diluted sample on Insulin ELISA (Exocell). Calculate ng/ml with appropriate dilution factor.
- 10) Run 20 μl of the neutralized solution in a Bradford Assay (250 μl Coomassie Blue Reagent, Thermo Scientific) against a standard curve. Calculate $\mu\text{g/ml}$ with appropriate dilution factor.
- 11) Divide Insulin content ng/ml by Protein content $\mu\text{g/ml}$.